

TTC7A Antibody



PACO12974

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TTC7A

Uniprot

Q9ULT0

Synonyms:

tetratricopeptide repeat domain 7A; TTC7A; MGC131720; MGC134830; TTC7 ;

Immunogen:

Human TTC7A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A