

TUBA1A Antibody



PACO12981

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TUBA1A

Uniprot

Q71U36

Synonyms:

50kd;

Immunogen:

Human TUBA1A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.320° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A