

TUBA8 Antibody



PACO12984

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TUBA8

Uniprot

Q9NY65

Synonyms:

449aa,50kd;

Immunogen:

Human TUBA8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.320° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A