

TUBB2B Antibody



PACO12986

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TUBB2B

Uniprot

Q9BVA1

Synonyms:

50kd;

Immunogen:

Human TUBB2B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.320° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A