

TULP2 Antibody



PACO12991

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TULP2

Uniprot

O00295

Synonyms:

tubby like protein 2; TULP2; TUBL2 ;

Immunogen:

Human TULP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A