

TXNL1 Antibody



PACO13003

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TXNL1

Uniprot

O43396

Synonyms:

thioredoxin-like 1; TXNL1; TRP32; TXL-1; TXNL; Txl ;

Immunogen:

Human TXNL1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A