

TYW1B Antibody



PACO13008

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TYW1B

Uniprot

Q6NUM6

Synonyms:

tRNA-yW synthesizing protein 1 homolog B (non-protein coding); TYW1B; MGC87315; NCRNA00069; RSAFD2 ;

Immunogen:

Human TYW1B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A