

UBASH3B Antibody



PACO13017

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UBASH3B

Uniprot

Q8TF42

Synonyms:

ubiquitin associated and SH3 domain containing, B; UBASH3B; KIAA1959; MGC15437; STS-1; STS1; p70 ;

Immunogen:

Human UBASH3B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A