

UBE2E3 Antibody



PACO13027

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

UBE2E3

Uniprot

Q969T4

Synonyms:

ubiquitin-conjugating enzyme E2E 3 (UBC4/5 homolog, yeast); UBE2E3; UBCH9; UbcM2 ;

Immunogen:

Human UBE2E3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A