

UBE2F Antibody



PACO13028

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UBE2F

Uniprot

Q969M7

Synonyms:

ubiquitin-conjugating enzyme E2F (putative); UBE2F; MGC18120; NCE2 ;

Immunogen:

Human UBE2F.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A