

UBE2G1 Antibody



PACO13029

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UBE2G1

Uniprot

P62253

Synonyms:

ubiquitin-conjugating enzyme E2G 1 (UBC7 homolog, yeast); UBE2G1; E217K; UBC7; UBE2G ;

Immunogen:

Human UBE2G1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A