

UBE2L6 Antibody



PACO13034

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

UBE2L6

Uniprot

O14933

Synonyms:

ubiquitin-conjugating enzyme E2L 6; UBE2L6; MGC40331; RIG-B; UBCH8 ;

Immunogen:

Human UBE2L6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A