

UBE2Q2 Antibody



PACO13037

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UBE2Q2

Uniprot

Q8WVN8

Synonyms:

ubiquitin-conjugating enzyme E2Q family member 2; UBE2Q2; DKFZp762C143 ;

Immunogen:

Human UBE2Q2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A