

UBE2V1 Antibody



PACO13041

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UBE2V1

Uniprot

Q13404

Synonyms:

ubiquitin-conjugating enzyme E2 variant 1; UBE2V1; CIR1; CROC-1; CROC1; UBE2V; UEV-1; UEV1; UEV1A ;

Immunogen:

Human UBE2V1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A