

# UBTD2 Antibody



PACO13058

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

UBTD2

**Uniprot**

Q8WUN7

**Synonyms:**

ubiquitin domain containing 2; UBTD2; DC-UbP; DCUBP; MGC30022 ;

**Immunogen:**

Human UBTD2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A