

UBXN11 Antibody



PACO13060

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UBXN11

Uniprot

Q5T124

Synonyms:

UBX domain protein 11; UBXN11; COA-1; DKFZp686F04228; PP2243; SOC; SOCI; UBXD5 ;

Immunogen:

Human UBXN11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A