

UCK1 Antibody



PACO13066

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UCK1

Uniprot

Q9HA47

Synonyms:

uridine-cytidine kinase 1; UCK1; FLJ12255; URK1 ;

Immunogen:

Human UCK1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A