

UFSP2 Antibody



PACO13074

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UFSP2

Uniprot

Q9NUQ7

Synonyms:

UFM1-specific peptidase 2; UFSP2; C4orf20; FLJ11200 ;

Immunogen:

Human UFSP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A