

UGDH Antibody



PACO13078

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UGDH

Uniprot

O60701

Synonyms:

UDP-glucose dehydrogenase; UGDH; GDH; UDP-GlcDH; UDPGDH; UGD ;

Immunogen:

Human UGDH.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A