

UGP2 Antibody



PACO13079

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UGP2

Uniprot

Q16851

Synonyms:

UDP-glucose pyrophosphorylase 2; UGP2; UDPG; UDPGP2; UGPP2; pHC379 ;

Immunogen:

Human UGP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A