

UGT3A1 Antibody



PACO13083

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UGT3A1

Uniprot

Q6NUS8

Synonyms:

UDP glycosyltransferase 3 family, polypeptide A1; UGT3A1; FLJ34658 ;

Immunogen:

Human UGT3A1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A