

UHMK1 Antibody



PACO13086

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UHMK1

Uniprot

Q8TAS1

Synonyms:

U2AF homology motif (UHM) kinase 1; UHMK1; KIS; Kist ;

Immunogen:

Human UHMK1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A