

UNC45B Antibody



PACO13095

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UNC45B

Uniprot

Q8IWX7

Synonyms:

unc-45 homolog B (C. elegans); UNC45B; CMYA4; FLJ38610; MGC119540; MGC119541; SMUNC45; UNC45 ;

Immunogen:

Human UNC45B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A