

UNG Antibody



PACO13096

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UNG

Uniprot

P13051

Synonyms:

uracil-DNA glycosylase; UNG; DGU; DKFZp781L1143; HIGM4; UDG; UNG1; UNG15; UNG2 ;

Immunogen:

Human UNG.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A