

# UNKL Antibody



PACO13097

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

UNKL

**Uniprot**

Q9H9P5

**Synonyms:**

unkempt homolog (Drosophila)-like; UNKL; C16orf28; FLJ12623; FLJ23360; KIAA0734; MGC5179; ZC3H5L; ZC3HDC5L ;

**Immunogen:**

Human UNKL.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A