

UQCR11 Antibody



PACO13100

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UQCR11

Uniprot

O14957

Synonyms:

ubiquinol-cytochrome c reductase, 6.4kDa subunit; UQCR;0710008D09Rik; QCR10 ;

Immunogen:

Human UQCR.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A