

UQCRQ Antibody



PACO13106

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UQCRQ

Uniprot

O14949

Synonyms:

ubiquinol-cytochrome c reductase, complex III subunit VII, 9.5kDa; UQCRQ; QCR8; QP-C; QPC ;

Immunogen:

Human UQCRQ.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A