

# URM1 Antibody



PACO13107

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

URM1

**Uniprot**

Q9BTM9

**Synonyms:**

ubiquitin related modifier 1 homolog (S. cerevisiae); URM1; C9orf74; MGC2668 ;

**Immunogen:**

Human URM1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A