

PACO13112

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

USMG5

Uniprot

Q96IX5

Synonyms:

up-regulated during skeletal muscle growth 5 homolog (mouse); USMG5; DKFZp566D211; HCVFTP2; MGC14697; bA792D24.4 ;

Immunogen:

Human USMG5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A