

USP20 Antibody



PACO13123

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

USP20

Uniprot

Q9Y2K6

Synonyms:

ubiquitin specific peptidase 20; USP20; KIAA1003; LSFR3A; VDU2 ;

Immunogen:

Human USP20.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A