

UTP14A Antibody



PACO13145

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UTP14A

Uniprot

Q9BVJ6

Synonyms:

UTP14, U3 small nucleolar ribonucleoprotein, homolog A (yeast); UTP14A; KIAA0266; NY-CO-16; SDCCAG16; dJ537K23.3 ;

Immunogen:

Human UTP14A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A