

UTP15 Antibody



PACO13146

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UTP15

Uniprot

Q8TED0

Synonyms:

UTP15, U3 small nucleolar ribonucleoprotein, homolog (S. cerevisiae); UTP15; FLJ12787 ;

Immunogen:

Human UTP15.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A