

UTP3 Antibody



PACO13149

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

UTP3

Uniprot

Q9NQZ2

Synonyms:

UTP3, small subunit (SSU) processome component, homolog (S. cerevisiae); UTP3; CRL1; CRLZ1; FLJ23256; SAS10 ;

Immunogen:

Human UTP3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A