

VAMP1 Antibody



PACO13154

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VAMP1

Uniprot

P23763

Synonyms:

vesicle-associated membrane protein 1 (synaptobrevin 1); VAMP1; DKFZp686H12131; SYB1; VAMP-1 ;

Immunogen:

Human VAMP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A