

# VAPA Antibody



PACO13160

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

VAPA

**Uniprot**

Q9P0L0

**Synonyms:**

VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa; VAPA; MGC3745; VAP-33; VAP-A; VAP33; hVAP-33 ;

**Immunogen:**

Human VAPA.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A