

VAV1 Antibody



PACO13166

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

VAV1

Uniprot

P15498

Synonyms:

vav 1 guanine nucleotide exchange factor; VAV1; VAV ;

Immunogen:

Human VAV1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A