

# VCPIP1 Antibody



PACO13170

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

VCPIP1

**Uniprot**

Q96JH7

**Synonyms:**

valosin containing protein (p97)/p47 complex interacting protein 1; VCPIP1; DKFZp686G038; DUBA3; FLJ23132; FLJ60694; KIAA1850; VCIP135 ;

**Immunogen:**

Human VCPIP1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A