

VDAC2 Antibody



PACO13172

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

VDAC2

Uniprot

P45880

Synonyms:

voltage-dependent anion channel 2; VDAC2; FLJ23841 ;

Immunogen:

Human VDAC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A