

VGLL1 Antibody



PACO13176

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

VGLL1

Uniprot

Q99990

Synonyms:

vestigial like 1 (Drosophila); VGLL1; TDU; VGL1 ;

Immunogen:

Human VGLL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A