

VMA21 Antibody



PACO13179

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VMA21

Uniprot

Q3ZAQ7

Synonyms:

VMA21 vacuolar H⁺-ATPase homolog (*S. cerevisiae*); VMA21; MGC125514; MGC125516; MGC131652 ;

Immunogen:

Human VMA21.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A