

VN1R4 Antibody



PACO13181

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VN1R4

Uniprot

Q7Z5H5

Synonyms:

vomeronasal 1 receptor 4; VN1R4; V1RL4 ;

Immunogen:

Human VN1R4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A