

VNN1 Antibody



PACO13182

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VNN1

Uniprot

O95497

Synonyms:

vanin 1; VNN1; HDLCQ8; MGC116930; MGC116931; MGC116932; MGC116933; Tiff66 ;

Immunogen:

Human VNN1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A