

VPREB3 Antibody



PACO13184

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VPREB3

Uniprot

Q9UKI3

Synonyms:

pre-B lymphocyte 3; VPREB3;8HS20; N27C7-2 ;

Immunogen:

Human VPREB3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A