

VPS33A Antibody



PACO13192

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

VPS33A

Uniprot

Q96AX1

Synonyms:

vacuolar protein sorting 33 homolog A (S. cerevisiae); VPS33A; FLJ22395; FLJ23187 ;

Immunogen:

Human VPS33A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A