

VPS33B Antibody



PACO13193

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

VPS33B

Uniprot

Q9H267

Synonyms:

vacuolar protein sorting 33 homolog B (yeast); VPS33B; FLJ14848 ;

Immunogen:

Human VPS33B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A