

VPS52 Antibody



PACO13198

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

VPS52

Uniprot

Q8N1B4

Synonyms:

vacuolar protein sorting 52 homolog (S. cerevisiae); VPS52; ARE1; DKFZp5471194; RP5-1033B10; SAC2; SACM2L; dJ1033B10.5 ;

Immunogen:

Human VPS52.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A