

VPS8 Antibody



PACO13202

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

VPS8

Uniprot

Q8N3P4

Synonyms:

vacuolar protein sorting 8 homolog (S. cerevisiae); VPS8; FLJ32099; KIAA0804 ;

Immunogen:

Human VPS8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A