

VSIG1 Antibody



PACO13205

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

VSIG1

Uniprot

Q86XK7

Synonyms:

V-set and immunoglobulin domain containing 1; VSIG1;1700062D20Rik; GPA34; MGC44287; dJ889N15.1 ;

Immunogen:

Human VSIG1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A