

# VTI1A Antibody



PACO13208

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

VTI1A

**Uniprot**

Q96AJ9

**Synonyms:**

vesicle transport through interaction with t-SNAREs homolog 1A (yeast); VTI1A; MVti1; Vti1-rp2;

**Immunogen:**

Human VTI1A.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A