

WBSCR16 Antibody



PACO13221

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RCC1L

Uniprot

Q96I51

Synonyms:

Williams-Beuren syndrome chromosome region 16; WBSCR16; DKFZp434D0421; MGC189739; MGC44931 ;

Immunogen:

Human WBSCR16.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A